# MEDsan® SARS-CoV-2 Antigen Rapid Test

Please read this package insert carefully prior to use and strictly follow the instructions.

## INTENDED USE

MEDsan® SARS-CoV-2 Antigen Rapid Test is a solid phase immunochromatographic assay intended for the *in vitro* qualitative detection of specific severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) antigens in human nasopharyngeal and/or oropharyngeal secretion. The test kit is applicable in healthcare system and the scientific field of research and for professional use only. Antigen is generally detectable in nasopharyngeal or oropharyngeal secretion during the acute phase of infection. Positive results indicate the presence of viral antigens, but clinical correlation with patient history and other diagnostic information is necessary to determine infection status. Negative results do not rule out SARS-CoV-2 infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19 and confirmed with a molecular assay, if necessary.

## INTRODUCTION

Coronaviruses are enveloped RNA viruses that are distributed broadly among humans and cause respiratory, enteric, hepatic, and neurologic diseases. Seven coronavirus species are known to cause human diseases. Four of these viruses - 229E, OC43, NL63, and HKU1 - are prevalent and typically cause common cold symptoms in immunocompetent individuals. The three other strains - severe acute respiratory syndrome coronavirus (SARS-CoV), Middle East respiratory syndrome coronavirus (MERS-CoV) and severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) - are zoonotic in origin and have been linked to sometimes fatal illness. Coronavirus disease 2019 (COVID-19) is a respiratory infectious disease caused by SARS-CoV-2. The most common symptoms include fever, cough, fatigue, shortness of breath, and loss of smell and taste. Nasal congestion, runny nose, sore throat, myalgia and diarrhea are found in a few cases. Currently, persons infected by SARS-CoV-2 are the main source of transmission. Asymptomatic infected people can also spread the virus. Based on the current epidemiological investigation, the incubation period is 2 to 14 days (median incubation time around 5 days).

## PRINCIPLE

The MEDsan® SARS-CoV-2 Antigen Rapid Test is a lateral flow immunochromatographic assay. The test uses SARS-CoV-2 antibodies (test line T) and goat anti-mouse IgG (control line C) immobilized on a nitrocellulose strip. The burgundy colored conjugate pad contains colloidal gold conjugated to SARS-CoV-2 antibodies (SARS-CoV-2 conjugates) and mouse IgG-gold conjugates. When a specimen followed by assay diluent is added to the sample well, SARS-CoV-2 antigen, if present, will bind to SARS-CoV-2 conjugates forming antigen antibodies complex. This complex migrates through nitrocellulose membrane by capillary action. When the complex meets the line of the corresponding immobilized antibodies, the complex will be combined forming a burgundy colored band which confirms a reactive test result. Absence of a colored band in the test region indicates a non-reactive test result. In addition, the test contains an internal control (C band) which should exhibit a burgundy colored band of the immunocomplex goat anti-mouse IgG/mouse IgG-gold conjugate regardless of the color development on any of the test bands. Otherwise, the test result is invalid and the specimen must be retested with another device.

## MATERIALS SUPPLIED

Each sealed pouch contains a test device and a desiccant. Buffer(s) | Sterile and single use specimen collection swab(s) | Single use extraction tube(s) with integrated dispensing tip(s) | Package insert

## MATERIAL REQUIRED BUT NOT PROVIDED

Materials not supplied but recommended for the performance are personal protection, as gloves and mouth protection. Standard microbiological supplies and equipment such as timer and tongue depressor (only required for oropharyngeal specimen collection) are not provided. External positive and negative controls can be purchased separately from MEDsan GmbH. These should be tested periodically consistent with good laboratory practice.

## STORAGE AND STABILITY

The test kit should be stored in a dry place protected from direct sunlight at 2-30 °C. The test device should be used within 1 hour after opening of the sealed pouch. If in a high humidity environment, use it immediately. DO NOT FREEZE. Do not use after the expiry date.

## WARNINGS AND PRECAUTIONS

- For *in vitro* diagnostic use only.
- The test is for single use only. Do not reuse.
- Do not perform the test in a room with strong air flow and in environment that is too hot, too humid, or too dry.
- The test device should be used as soon as possible after opening the pouch. Avoid keeping it in the air for a long time, which may result in failure due to damp. Do not use it if the pouch is damaged or broken.
- This test is only validated using the material provided with this kit.
- Do not mix components from different lots.
- Handle all specimens as if infectious by using safe laboratory procedures.
- When testing many samples of specimens, please mark well to avoid confusion.
- After the test is completed, used materials as test device, extraction tube and swabs should be discarded into medical waste garbage bags, which will be specially disposed by the qualified unit to handle medical waste.
- This test has been authorized only for the detection of SARS-CoV-2 proteins, not for any other viruses or pathogens.

## SPECIMEN COLLECTION

Standard precautions should always be followed whenever samples are obtained from patients: use protective gown, pair of nonsterile gloves, face mask and visor for face and eye protection. Prepare the extraction tube (refer to section 'Test Procedure') and use the supplied sterile, single use specimen collection swabs.

### Nasopharyngeal specimen collection

- Ask the patient to take off the mask and to blow their nose to clear nasal passage of excessive mucus.
- Tilt the patient's head back 70 degrees.
- Insert swab into the nostril. Swab should reach the surface of posterior nasopharynx. Gently swirl the swab for 5-10 seconds to absorb secretions.
- Gently remove swab while rotating it.
- Place the swab into prepared extraction tube (refer to point 2 in section 'Test Procedure').
- Ask the patient to reapply the mask.

### Oropharyngeal specimen collection

- Ask the patient to take off the mask.
- Tilt the patient's head back 70 degrees.
- Insert swab into the nostril. Swab should reach the surface of posterior nasopharynx. Gently swirl the swab for 5-10 seconds to absorb secretions.
- Place the swab into prepared extraction tube (refer to point 2 in section 'Test Procedure').
- Ask the patient to reapply the mask.

### Combination of oropharyngeal and nasopharyngeal specimen collection

It is also possible to combine both collection methods using the same swab (first oropharyngeal, then nasopharyngeal).

English


MEDsan®  
biological health solutions

## TEST PROCEDURE

This test procedure has to be read completely before performing the test.

Allow test device, specimen, buffer and/or controls to equilibrate to room temperature (15-30 °C) prior to testing and use it as soon as possible.

- Add all of one single buffer (380 µL) into a single use extraction tube.
- After specimen collection (refer to section 'Specimen Collection') insert the swab into the extraction tube which contains 380 µL of the buffer and rotate the swab constantly. Repeat several times and incubate for at least 1 minute.
- Squeeze the swab on the tube wall so that the liquid is screwed out. Take out and discard the swab according to the treatment of medical waste.
- Cover the extraction tube with its dispensing tip.
- Remove the test device from the sealed foil pouch and place it on a clean and even surface. Add 2 drops of the sample solution vertically into the sample well of the test device.
- Wait 15 minutes to interpret and record the test result. The result is invalid after 20 minutes.

## INTERPRETATION OF RESULTS

### POSITIVE


If the C band and T band are both present, then the test indicates the presence of SARS-CoV-2 antigens in the specimen. The test result is positive. The purple red test line may vary in shade and intensity depending on the detected antigen concentration. Also, a light or faint test line must be interpreted as a positive result.

### NEGATIVE

If only the C band is present, the absence of any burgundy color in the T band indicates that no SARS-CoV-2 antigens are detected in the specimen. The test result is negative.

### INVALID

Control line C is missing, or control line C and test line T are missing. Incorrect specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Read the instructions carefully again and repeat the test with a new test device. If the problem persists, discontinue using the test device immediately and contact your local distributor.

## LIMITATIONS

- This test is suitable for testing human nasopharyngeal and/or oropharyngeal secretion. This test kit is not intended to be used for other body fluids and samples.
- The test results should be used in combination with the clinical examination, medical history, and other examination results.
- A negative result for an individual subject indicates absence of detectable SARS-CoV-2 antigens. A negative test result does not preclude the possibility of exposure to or infection with SARS-CoV-2.
- A negative result may occur if the quantity of the SARS-CoV-2 antigens present in the specimen is below the detection limits of the assay.
- Positive test results do not rule out co-infections with other pathogens.
- Positive test results do not differentiate between SARS-CoV and SARS-CoV-2.
- Negative test results are not intended to rule out other non-SARS viral or bacterial infections.
- Optimal test performance requires strict compliance with the test procedure described in this instructions for use. Deviations may lead to aberrant results. Incorrect specimen volume may lead to invalid test results.
- Do not keep your prepared sample solution longer than for 60 minutes. This may lead to false test results.

## PERFORMANCE CHARACTERISTICS

### 1. Clinical Studies

Reagents have been sent to independent laboratories for clinical evaluation. Antigen detection in the samples of COVID-19 patients has a high consistency with nucleic acid detection from swab samples. The relative sensitivity is 92.5% (96.5% for samples with Ct values ≤ 33). The relative specificity rate is 99.8% and the accuracy of the product is 98.4%.

### 2. Limit of Detection (LoD)

Limit of detection (LoD) of the MEDsan® SARS-CoV-2 Antigen Rapid Test was determined by two different methods. During the first method different concentrations of heat inactivated SARS-CoV-2 were evaluated. The MEDsan® SARS-CoV-2 Antigen Rapid Test is confirmed with a LoD of 14.4 TCID<sub>50</sub>/mL. The second method used different concentrations of recombinant antigens demonstrating a LoD of 10 pg/mL.

### 3. Analytical Specificity/Cross Reactivity/Microbial Interference

Analytical specificity of the MEDsan® SARS-CoV-2 Antigen Rapid Test has been evaluated to other pathogens. No antigen false positive results or microbial interferences were observed with the following potential cross reactants: Human Coronavirus (229E, OC43, NL63, HKU1), MERS-CoV, Adenovirus, Human Metapneumovirus (hMPV), Parainfluenza virus 1-4, Influenza A, Influenza B, Influenza C, Enterovirus, Respiratory syncytial virus, Rhinovirus, *Haemophilus influenzae*, *Streptococcus pneumoniae*, *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Streptococcus pyogenes*, *Candida albicans*, *Bordetella pertussis*, *Mycoplasma pneumoniae*, *Chlamydia pneumoniae*, *Legionella pneumophila*, and a pooled human nasal wash – representative of normal respiratory microbial flora (healthy donors).

### 4. Interference Substances Studies

Potential interference of the MEDsan® SARS-CoV-2 Antigen Rapid Test was evaluated using natural clinical samples. No antigen false negative or false positive results have been observed with the following potential interference substances at the stated concentrations: human blood (1% v/v), mucosal protein (1 mg/mL), menthol (50 mg/mL), dyclonine/menthol (2 mg/mL), phenylephrine (1% v/v), oxymetazoline (1% v/v), triamcinolone (50 mg/mL), ribavirin (50 mg/L), alkalol (10% v/v), benzocaine and menthol (50 mg/mL), fluticasone propionate (5% v/v), tobramycin (8 µg/mL), mupirocin (10 mg/mL), and biotin (0.15 mg/mL).

## REFERENCES

- Weiss SR, Leibowitz JL. Coronavirus pathogenesis. Adv Virus Res 2011; 81: 85-164.
- Masters PS, Perlman S. Coronaviidae. In: Knipe DM, Howley PM, eds. Fields virology. 6th ed. Lippincott Williams & Wilkins, 2013: 825-88.
- Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. Trends Microbiol 2016; 24: 490-502.
- Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. Nat Rev Microbiol 2019; 17: 181-192.
- World Health Organization. Naming the Coronavirus disease (COVID-19) and the virus that causes it. World Health Organization. Archived from the original on 28 February 2020. Retrieved 28 February 2020.
- Hessen MT (27 January 2020). "Novel Coronavirus Information Center: Expert guidance and commentary". Elsevier Connect. Archived from the original on 30 January 2020. Retrieved 31 January 2020.
- CDC Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens for COVID-19. <https://www.cdc.gov/coronavirus/2019-nCoV/lab/guidelines-clinical-specimens.html>.

## INDEX OF SYMBOLS

	See instruction for use		Expiry date
	For <i>in vitro</i> diagnostic use only		Batch number
	Store between 2-30 °C		Manufacturer
	Tests per kit		Keep dry
	Catalog number		Do not reuse
	Keep away from sunlight		

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REFERENZEN

1. Weiss SR, Leibowitz JL. Coronavirus pathogenesis. Adv Virus Res 2011; 81: 85-164.
2. Masters PS, Perlman S. Coronaviridae. In: Knipe DM, Howley PM, eds. Fields virology. 6th ed. Lippincott Williams & Wilkins, 2013: 825-58.
3. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. Trends Microbiol 2016; 24: 490-502.
4. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. Nat Rev Microbiol 2019; 17: 181-192.
5. „Naming the Coronavirus disease (COVID-19) and the virus that causes it“. World Health Organization. Archived from the original on 28 February 2020. Retrieved 28 February 2020.
6. Hessen MT (27 January 2020). „Novel Coronavirus Information Center: Expert guidance and commentary“. Elsevier Connect. Archived from the original on 30 January 2020. Retrieved 31 January 2020.
7. CDC Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens for COVID-19 https://www.cdc.gov/coronavirus/2019-nCoV/lab/guidelines-clinical-specimens.html

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REFERENSER

1. Weiss SR, Leibowitz JL. Coronavirus pathogenesis. Adv Virus Res 2011; 81: 85–164.
2. Masters PS, Perlman S. Coronaviridae. In: Knipe DM, Howley PM, eds. Fields virology, 6th ed. Lippincott Williams & Wilkins, 2013: 825–58.
3. Su S, Wong G, Shi W, et al. Epidemiology, genetic recombination, and pathogenesis of coronaviruses. Trends Microbiol 2016; 24: 490–502.
4. Cui J, Li F, Shi ZL. Origin and evolution of pathogenic coronaviruses. Nat Rev Microbiol 2019; 17: 181–192.
5. „Naming the Coronavirus disease (COVID-19) and the virus that causes it“. Världshälsoorganisationen. Arkiverad från originalet den 28 februari 2020. Hämtad den 28 februari 2020.
6. Hsien MT (27 januari 2020). „Novel Coronavirus Information Center: Expert guidance and commentary“. Elsevier Connect. Arkiverad från originalet den 30 januari 2020. Hämtad den 31 januari 2020.
7. CDC Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens for COVID-19. https://www.cdc.gov/coronavirus/2019-nCoV/lab/guidelines-clinical-specimens.html

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